# **ROTORS, TUBES, BOTTLES AND CAPS**

FOR Hitachi Preparative Ultracentrifuge

This manual contains instruction for operation of Angle rotor, Neo-angle rotor, Swing rotor and Vertical rotor. Carefully read this manual in conjunction with manual for centrifuge and use the rotor correctly. Retain it for future reference after reading.

The names of actual companies and products mentioned herein may be the trademarks of their respective owners.

## ABOUT ⚠ MARKS

In this manual,  $\bigwedge$  marks are used to make attention, to prevent personal injury or machine damage. The meaning of the mark is as follows.

 $\underline{ \boldsymbol{\Lambda}} \ \ \text{WARNING}: \ \text{Indicates a potentially hazardous situation which, if not avoided, could result in severe }$ 

personal injury or possible death.

⚠ CAUTION: Indicates a hazardous situation which, if not avoided, could result in personal injury, or

severe damage to the instrument.

#### SAFETY REMINDER

Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly. For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

## **⚠** WARNING

- ◆ Never use any sample capable of producing flammable or explosive vapors. Your ultracentrifuge and rotor have no explosion-proof construction.
- ◆ Your ultracentrifuge and rotor are not designed to confine any sample particles dispersed due to a leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, take additional precautions to prevent exposure to these materials, (e.g., use of isolated areas).
- ◆ Never exceed the maximum speed of the rotor stamped on it: Exceeding the specified speed could break the rotor, resulting in damage to the ultracentrifuge.
- ◆ Check the chemical resistance chart attached the rotor, and do not use any sample inapplicable to the rotor (including buckets), Using such a sample could corrode the rotor (including buckets).
- ◆ Do not remove the RLM adapter or optical adapter from rotor, or replace it with the adapter of another rotor. The adapter is a critical component that detects the over-speed of rotor: If an adapter that is not compatible with the rotor is attached, the rotor could break, resulting in damage to the ultracentrifuge.
- ◆ Do not allow the rotor temperature to rise above 100°C: This would cause the material to become brittle.
- ◆ Never autoclave the rotor or sterilize it in boiling water: The strength of rotor may be significantly lost.
- When using swing rotor, be sure to set all buckets of the same type whether or not samples are put in buckets: Failure to do so could not only cause the rotor to vibrate, but could result in the rotor being deformed and the buckets becoming detached, which is very dangerous. Never use buckets made by other companies, or any other type of bucket that is not exclusively made for the rotor even if it is made by Hitachi Koki.
- If the centrifuge, rotor, or an accessory is contaminated by samples that toxic or radioactive, or blood samples that pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- ◆ If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from a Hitachi Koki authorized sales/service representative.
- ◆ It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to a Hitachi Koki authorized sales/service representative.

# **⚠** CAUTION

- ◆ Check the chemical resistance chart in separate manual, and do not use any sample inapplicable to the tubes, the bottles, the tube adapters, or the caps. Using such a sample could deteriorate them.
- ◆ The allowable speed may be lower depending on the density of sample, use of salt-such as cesium chloride-or the combination of tubes and bottles, or tube adapters and caps. Refer to this instruction manual.
- Always balance the tubes/bottles and samples within the allowable imbalance of rotor, including the tube adapters and caps: Operation with imbalance exceeding the allowable range could damage the ultracentrifuge and rotor.
- ◆ Do not load only one tube or load tubes asymmetrically: Asymmetrical loading may cause imbalance operation and damage the centrifuge and rotor.
- ◆ Clean the inside of drive hole (crown hole) in the rotor and the surface of drive shaft (crown) of centrifuge once a month. If the drive hole or drive shaft is stained or if any foreign matter adheres, the rotor may have been improperly installed and could come off during operation.
- Do not use tubes/bottles that have exceeded their life expectancy. Failure to do so could result in damage of tubes/bottles and the rotor and the centrifuge.
  - The life expectancy of tubes/bottles depends on factors such as the characteristics of samples, speed of the rotor used, and temperature.
  - Always check for deterioration and damage (cracks, deformation, and so on) on tubes/bottles before using them. Do not use the tubes/ bottles if you find such a problem.
- ◆ Maintain and inspect the rotor each time after using it. If you fine any abnormality in the rotor, stop using it and contact a Hitachi Koki sales/service representative.

# **CONTENTS**

1.	Specification of rotors	1
1	.1 Preface	1
1	.2 Angle rotors	2
1	.3 Neo-angle rotors	3
1	.4 Swing rotors	4
1	.5 Vertical rotors	5
1	.6 Prevention of over-speed	6
1	.7 Relation between rotors and ultracentrifuges	7
1	.8 Calculation of separation characteristics	9
2.	Tubes, bottles and caps	11
2	2.1 The kind of tubes, bottles and caps ·····	11
2	2.2 Materials of tubes·····	13
2	2.3 Washing ·····	13
2	2.4 Sterilization ·····	14
2	2.5 Storage·····	15
2	2.6 Inspection ·····	15
2	2.7 Life	16
2	2.8 How to use tubes and bottles·····	17
	Tubes, 5mL and over, for angle rotors and vertical rotors	
	Tubes, 4mL or less, for angle rotors	21
	• Tubes (12PA/12PE/40PA/40PE tubes) that are applicable for S-Cap	23
	Bottles	27
	Seal tubes	28
	Thick-walled tubes for angle rotors	30
	Tubes for swing rotors	31
	Re-seal tubes and cone-top tubes	32
2	2.9 How to balance ·····	35
3.	How to use rotors	36
3	3.1 Allowable rotor speeds ·····	36
3	3.2 Sample limitation ·····	37
3	3.3 Adhesion of sample, etc	37
3	3.4 Angle rotors ·····	38
3	3.5 Swing rotors ·····	40
3	8.6 Neo-angle rotors and titanium Vertical rotors ······	45

4. Mainter	nance48
4.1 Anti-co	prrosion property
4.2 Mainte	enance of rotors48
4.3 Steriliz	ring rotors 50
4.4 Inspec	ction of rotors50
4.5 Repla	cing optical adapter/disk51
4.6 Life of	rotors54
5. Trouble	eshooting 55
5.1 Rotors	5555
5.2 Tubes	56
Appendix 1	The list of caps for tubes
Appendix 2	The list of S-Series Cap and tools 60
Appendix 3	The list of space caps 62
Appendix 4	Tightening tools of the cap63
Appendix 5	The list of bottles
Appendix 6	The list of tube kits 65
Warranty	67
Decontamina	ation69
Rotor retiren	nent69

# 1. Specification of rotors

#### 1.1 Preface

The rotors can be classified into following matters:

#### (1) Classified into shapes

Angle rotors......This rotor holds the tubes and bottles at an angle to the axis of rotation. It is mainly used to separate components in a cell using fraction centrifugation (Pelleting) and separate nucleic acid using density gradient sedimentation equilibrium.

Neo-angle rotors.......This rotor holds the tubes at a smaller angle to the axis of rotation than the angle rotors. It is effective to separate samples that produce precipitates and floating matter such as plasmid DNA in a short time.

Vertical rotors .......This rotor holds the tubes perpendicularly to the centrifugal force. It is expected for separation in a shorter time with this rotor than with swing rotors and angle rotors. It is specially effective to separate nucleic acid using the density gradient sedimentation equilibrium.

#### (2) Classified into materials

Titanium rotors......Rotors made from titanium alloy.

Aluminum rotors .....Rotors made from aluminum alloy.

## (3) Classified into management type

Rotor with optical adapter/disk

.....The rotor's operation history is recorded in the log book for control. An optical adapter/disk is attached as an optical over-speed prevention device.

RLM rotor ......This is exclusively used with a CP-NX, a CP-WX, a CP-MX, or a CP-α ultra-centrifuge. The rotor life management (RLM) system controls the rotor's operation history automatically. An RLM adapter is attached for magnetic recording.

# 1.2 Angle rotors

Angle rotors are made from either aluminum alloy or titanium alloy. Table 1-1 shows the specifications of angle rotors and Fig. 1-1 shows the name of parts.




Fig. 1-1 Names of parts of angle rotor

Table 1-1 Specification of angle rotor

Rotor	Maximum Speed (rpm)	Maximum R.C.F (× g)	K factor	Tube* (mL× Number)	Rotor* capacity (mL)	Weight (kg)	Rotor body material	Cover material
P100AT2	100,000	803,000	18	6.5 × 8	52	3.7	Titanium alloy	Titanium alloy
P90AT	90,000	700,000	25	12 x 8	96	4.7	Titanium alloy	Titanium alloy
P80AT	80,000	615,000	27	12× 8	96	6.0	Titanium alloy	Titanium alloy
P70AT2	70,000	452,000	36	12 × 12	144	6.7	Titanium alloy	Titanium alloy
P70AT	70,000	505,000	44	40 × 8	320	10.0	Titanium alloy	Titanium alloy
P65A	65,000	370,000	48	12 × 10	120	4.1	Aluminum alloy	Aluminum alloy
P50AT4	50,000	316,000	32	6.5 × 44	286	11.1	Titanium alloy	Aluminum alloy
P50AT2	50,000	303,000	70	40 × 12	480	15.6	Titanium alloy	Aluminum alloy
P50A3	50,000	252,000	33	1.5 × 24	36	2.1	Aluminum alloy	Aluminum alloy
P45AT	50,000	235,000	130	94× 6	564	14.6	Titanium alloy	Aluminum alloy
P42AT	42,000	223,000	12	$0.23 \times 72$	16.6	8.0	Titanium alloy	Aluminum alloy
P32AT	32,000	111,000	186	12 × 32	384	8.8	Titanium alloy	Aluminum alloy
P27A	27,000	106,000	352	160 × 6	960	12.5	Aluminum alloy	Aluminum alloy
P21A2	21,000	71,000	486	230 × 6	1380	11.4	Aluminum alloy	Aluminum alloy
P19A	30,000	55,100	754	230 × 6	1380	18.7	Aluminum alloy	Aluminum alloy

<sup>\*</sup> Nominal capacity

# 1.3 Neo-angle rotors

Neo-angle rotors are made from titanium alloy. Table 1-2 shows the specifications of neo-angle rotors and Fig. 1-2 shows the name of parts.




Fig. 1-2 Names of parts of Neo-angle Rotor

Table 1-2 Specification of neo-angle rotor


	Rotor	Maximum Speed (rpm)	Maximum R.C.F (× g)	K factor	Tube* (mL×Number)	Rotor* capacity (mL)	Weight (kg)	Rotor body material	Cap material
	P90NT	90,000	646,000	10	5× 8	40	4.9	Titanium alloy	Aluminum alloy
	P65NT	65,000	402,000	23	12 × 10	120	9.6	Titanium alloy	Aluminum alloy
Ī	P65NT2	65,000	431,000	15	5 × 18	90	8.7	Titanium alloy	Aluminum alloy

<sup>\*</sup> Nominal capacity

## 1.4 Swing rotors

There are three types of swing rotors depending on buckets types (hook type buckets, pin type buckets (only P65ST), and top loading type buckets (only P32ST)). Table 1-3 shows the specifications of swing rotors and Fig. 1-3 shows the name of parts.

Swing rotors (hook type buckets and pin type buckets) other than the P32ST rotor



<sup>\*</sup> Some rotors do not have a rotor pin and a balance weight.

An old rotor may have a pin guide in place of a balance weight.

#### ●P32ST rotor (top loading type buckets)




Fig. 1-3 Names of parts of swing rotor

Table 1-3 Specification of swing rotor

Rotor	Maximum Speed (rpm)	Maximum R.C.F (x g)	K factor	Tube* (mL × Number)	Rotor* capacity (mL)	Weight (kg)	Rotor body material	Bucket material
P65ST	65,000	419,000	48	5×3	15	4.3	Titanium alloy	Titanium alloy
P56ST	56,000	409,000	54	4×6	24	5.2	Titanium alloy	Titanium alloy
P55ST2	55,000	366,000	50	5×6	30	5.5	Titanium alloy	Titanium alloy
P40ST	40,000	284,000	139	13×6	78	6.9	Titanium alloy	Titanium alloy
P32ST	32,000	180,000	198	40×6	240	7.1	Titanium alloy	Titanium alloy
P28S	28,000	141,000	252	40 × 6	240	6.2	Aluminum alloy	Titanium alloy

<sup>\*</sup> Nominal capacity

# 1.5 Vertical rotors

Vertical rotors are made from titanium alloy. Table 1-4 shows the specifications of vertical rotors and Fig. 1-4 shows the name of parts.




Fig. 1-4 Names of parts of vertical rotor

Table 1-4 Specification of vertical rotor

Rotor	Maximum Speed (rpm)	Maximum R.C.F (x g)	K factor	Tube* (mL × Number)	Rotor* capacity (mL)	Weight (kg)	Rotor body material	Cap or adapter material
P100VT	100,000	700,000	6	5× 8	40	4.1	Titanium alloy	Aluminum alloy
P65VT3	65,000	402,000	13	12 × 10	120	10.1	Titanium alloy	Aluminum alloy
P65VT2	65,000	416,000	10	5×16	80	9.5	Titanium alloy	Aluminum alloy
P50VT2	50,000	243,000	36	40 × 8	320	11.9	Titanium alloy	Aluminum alloy

<sup>\*</sup> Nominal capacity

## 1.6 Prevention of over-speed

If the rotor speed exceeds the maximum limit, the rotor might break, leading to a serious accident. The RLM adapter or optical adapter/disk fitted at the bottom of the rotor prevents overspeed of the rotor.

MARNING: Do not remove the RLM adapter or optical adapter from rotor, or replace it with the adapter for another rotor.

> The adapter is a critical component that detects the over-speed of rotor: If an adapter that is not compatible with the rotor is attached, the rotor could break, resulting in damage to the ultracentrifuge.

- ⚠ CAUTION: •Do not bring the RLM adapter near a magnet, or scratch it: Doing so will erase the memory stored in adapter, and make the rotor unusable. To prevent the adapter from being scratched, store the rotor with RLM adapter, using the stand provided with the rotor (rotor stand for protecting adapter). If the adapter is scratched, immediately stop using the rotor, and contact a Hitachi Koki authorized sales/service representative.
  - Take great care not to scratch the optical adapter: Any scratch on adapter will make the rotor unusable. To prevent the adapter from being scratched, store the rotor with optical adapter, using the stand provided with the rotor (rotor stand for protecting adapter). If the adapter is scratched, immediately stop using the rotor, and contact a Hitachi Koki authorized sales/service representative.
  - •If foreign matter adheres to the RLM adapter or optical adapter/disk, wipe it off promptly. Otherwise, the RLM adapter or optical adapter/disk may be corroded.

The RLM adapter of RLM Rotor is a magnetic memory. Use care to protect it from damage, and keep it away from magnetism (otherwise memory contents would be destroyed). To protect the RLM adapter, be sure to keep the RLM Rotor on the rotor stand provided.

The optical adapter/disk of the standard rotor has stripes of black (not reflecting light) and white (reflecting light) as shown in Fig.1-5 and the number of stripes depends on the maximum speed of the rotor.

This adapter/disk keeps the rotor speed under the maximum permissible value.

 See | 4.5 Replacing Optical Adapter/Disk for replacing the optical disk.




Fig. 1-5 Optical adapter/disk

At model 55PA, 65P and 55P-2, if over speeding occurs, a small pin jumps out from the pin guide and the overspeed preventing system functions to reduce speed.

When using an old type ultracentrifuge, refer to

"1.7 Relation between Rotor and Ultracentrifuge".




Fig. 1-6

## 1.7 Relation between rotors and ultracentrifuges

The system to prevent the over speed is different from the system of the present centrifuge. Therefore, usable rotor is limited. Use the proper rotor in accordance with Table 1-5.

**CP NX Series** Type of CP100NX CP90NX CP80NX Centrifuge **CP WX Series** 55P-3,80P,7Series(80P-7 etc.), CP100WX 72Series(85P-72 etc.), CP90WX 55PA Rotor 55P SCP-H Series(SCP85H etc.), CP80WX 55P-2 CP MX Series SCP-H2 Series(SCP85H2 etc.), 40P 65P CP100MX CP Series(CP70 etc.), CP56G CP80MX CP β Series(CP85 β etc.) CP70MX CP α Series CP100α CP90H α CP80H α / Note 1) Note 5) X X X Δ **RLM** adapter Note 2) Note 1) Note 3) 0 X Δ Note 4)  $\triangle$ Optical adapter/disk

Table 1-5 The relation between Rotor and Ultracentrifuge

 $\bigcirc$  :Usable  $\triangle$ :Refer to the notice  $\mathbf{X}$ :Impossible to use

#### Note 1)

Do not use a P21A2 rotor with the centrifuge other than a series of CP-NXWX/MX ultracentrifuges.

#### Note 2)

Following rotors shall not be used with the former centrifuge, model 55PA, 55P- 2 or 65P, to use them safely.

P70AT(RP70T), P50AT4(RP50AT4), P50AT2(RP50T- 2), P45AT(RP45T), P65NT2
P65NT(RP65NT), RPV65T, P65VT2(RP65VT2), P65VT3(RP65VT3), P50VT(RPV50T), RPV45T, P50VT2

#### Note 3)

Following rotors shall not be used with the former centrifuge model 55PA because the rotor chamber of model 55PA is small.

P1 9A(RP19), P40ST(RPS40T), P28S(SRP28SA), P28S2(SRP28SA1)

#### Note 4)

All of following rotors have the balance weight. In case of using the these rotors with the former centrifuge, model 55PA, 55P- 2 or 65P, install the pin guide to the rotor in place of the balance weight. In case the pin guide is necessary, please buy it.

DO NOT use rotor without pin guide with these centrifuge because it causes danger.

#### Note5)

The P32ST rotor is a RLM rotor. The P32ST rotor is not a rotor with an optical adapter/disk. The P32ST rotor can be used only in a series of CP-NXWX/MX ultracentrifuges.

P55AT(RP55T), SRP50AT, RP50T, P50A2(RP50-2), P42AT(RPL42T), P42A(RP42), RP40-2, RP40, RP30-3, P30A2(RP30-2), RP21, P19A(RP19), SRP28SA, SRP28SA1, RPW50T, RPW45, RPW35, P48ZT(RPZ48T), P35ZT(RPZ35T), P32CT(RPC32T)

Pin guide Balance weight

Fig. 1-7 Pin guide and balance weight

- Old rotors with black adapter are used with model 55PA, 55P- 2 or 65P only.
- Old rotors which have no adapter are used with model 40 P or 55P only.




Fig. 1-8

## 1.8 Calculation of separation characteristics

#### (1) Relationship between the relative centrifugal force and speed

The relative centrifugal force (RCF) generated by rotation is related to the speed and radial distance from the axis of rotation. The RCF is obtained from the following formula.

(Generally, the RCF is expressed as a ratio of the earth's gravitational acceleration and " x g" is used as the unit.)

$$RCF = 1.118 \times 10^{-5} \times r \times N^{2} (\times g)$$

N: Speed (rpm)

r : Radial distance from the axis of rotation (cm)

The speed (rpm) can also be calculated from the RCF ( x g) using this formula.

#### (2) To find the separation time

The sedimentation time T required for the particles scattered in the solution to lower from the top (r<sub>t</sub>) to the bottom (r<sub>b</sub>) of the tube. The sedimentation characteristic is expressed as sedimentation coefficient s but generally the sedimentation coefficients becomes a very small value and  $s \times 10^{13}$  may be indicated as S. This S expressed as  $s \times 10^{13}$  is called the sedimentation coefficient in Svedberg units. However, in biochemistry, etc., it may simply be called the sedimentation coefficient.

$$T(hr) = \frac{10^{13}}{3600 \times s} \cdot \frac{\ln{(r_b)} - \ln{(r_t)}}{\omega^2}$$

ω: Angular velocity  $ω = \frac{2π}{60} \cdot N$ 

N: Speed(rpm)

rt : Distance from axis of rotation to top of solution

rb : Distance from axis of rotation to bottom of tube

Here, if it is defined that K = 
$$\frac{\ln{(r_b)} - \ln{(r_t)}}{\omega^2} \cdot \frac{10^{13}}{3600}$$
,

$$T(hr) = \frac{K}{S}$$

If S is unknown, the following formula (Stokes' law) is used for calculation.

$$S = \frac{d^2(\rho_2 - \rho_1)}{18\eta} \times 10^{13}$$

d : Diameter of lowering particle (cm)

ρ<sub>1</sub>: Density of solution surrounding particles (g/mL)

ρ<sub>2</sub>: Density of lowering particle (g/mL)

γ : Viscosity of solution surrounding particles (Poise)
 (The viscosity of water at 20°C is 0.01 Poise)

S in the pure water at 20°C is expressed as  $S_{20,w}$ .

K is called the K factor (clearing factor) which varies depending of the type of rotor and the speed. When the rotor is used at a certain speed, the K factor  $(K_N)$  can be obtained from the formula shown below.

$$K_N = K \left( \frac{N_{max}}{N} \right)^2$$

N : Speed used (rpm)

N<sub>max</sub>: Maximum speed (rpm)

 $K_N$ : K factor at the speed used

K: K factor at the maximum speed

# 2. Tubes, bottles and caps

## 2.1 The kind of tubes, bottles and caps

#### Tubes and bottles

#### ○Tubes (thin-walled tubes)

For angle rotors or vertical rotors, tubes must be filled up and used with caps. For swing rotors, tubes must filled to within 3mm from top of tubes and caps are not used. SS tubes can be used with optional volume.

#### Thick-walled tubes

Thick-walled tubes are used in angle rotors and swing rotors without caps. For angle rotors, these can be used with optional volume less than net volume. For swing rotors, fill tubes upto within 3mm from the top of tubes with a sample.

#### ○Seal tubes

Seal tubes are heat sealed and used with space caps in angle rotors, neo-angle rotors and vertical rotors. The tube sealer, model STF-1, model STF2, model STF3 and the tube rack are necessary to seal these tubes.

#### ○Re-seal tubes<sup>™</sup> and Cone-top tubes<sup>™</sup> (we discontinue to sell them.)

Re-seal tubes and cone-top tubes are sealed by screwing a plug into the tube and used in angle rotors, neo-angle rotors and vertical rotors with a plug and a crown. These tubes must be filled up the sample. Specified tube vise and tightening tools are necessary to seal these tubes.

\*Registered trademark of Seton Scientific Company U. S. A.




Fig. 2-1

#### ○ Bottles

A container of the screw cap type is called a bottle. Bottles with inner caps are for high-speed rotation. Bottles with Plastic (polyacetal) screw caps are classified as bottles (B) and those with metal (aluminum alloy) caps are classified as bottles (C). (See Fig.2-2 and appendix 5.)

When a 30PC bottle (C3), 30PC bottle (C), 30PA bottle (C), 70PC bottle (B), 70PC bottle (C), 70PA bottle (C), 80PC bottle (C2) and 80PC bottle (C) are used at more than 100,000 × g (max G-Force), the liquid should be over the shoulder of the bottle. Any volume which does not exceed the net volume is applicable at less than 100,000×g.




Fig. 2-2 Bottles


#### Caps and adapters

#### ○Caps

A Cap is used when a thin walled tube is used in an angle rotor or vertical rotor. Be careful that different caps have different assemblies (parts) and materials (see appendix 1). Tightening tools for caps are required (see appendix 4).

#### Space caps (for seal tubes)

Space caps prevent deformation of seal tubes during operation and are made from aluminum alloy. Be careful that different rotors use different space caps (see appendix 3).






Fig. 2-3 Caps

Crowns and plugs prevent deformation of re-seal tubes or cone-top tubes and leakage from tubes during operation. Crowns are made from aluminum alloy and plugs are made from stainless steel. Be different rotors use different crowns or plugs.

#### Adapters

Adapters are used when relatively smaller tubes are inserted into the tube cavities. Most adapters are made from white polyacetal. Since the specific gravity of adapters is high, the rotation speed is restricted when adapters are used.




Fig. 2-4 Adapters

#### 2.2 Materials of tubes

⚠ CAUTION: Run the rotor at the proper temperature (see the below). Otherwise the tubes may be damaged or deformed during operation. When running the rotor beyond this limit, test the tube to ensure that there is no abnormality before actual operation.

↑ CAUTION: Check the chemical resistance chart in separate manual, and do not use any sample inapplicable to the tubes and the bottles.

#### PC; Polycarbonate

Strong, transparent and can be autoclaved at 121°C or less. Use the tubes/bottles at a temperature between 4 °C and 25 °C. Weak in organic solution, alkali solution, and alkali detergent.

#### PET; Polyethylene terephthalate

Transparent. Can not be autoclaved. Sliceable and puncturable. Weak in organic solution. Use the tubes/bottles at a temperature between 4°C and 20°C.

#### PE; Polyethylene

Opaque. Can not be autoclaved. Excellent in chemical resistance.

Use the tubes/bottles at a temperature between 4°C and 20°C.

#### PP: Polypropylene

Strong but susceptible to low temperature (brittle point: 0°C). It can be autoclaved at 121°C or less.

Use the tubes/bottles at a temperature between 4°C and 25°C.

#### PA; Polypropylene copolymer

Strong and can be autoclaved at 115°C or less. Use the tubes/bottles at a temperature between 4°C and 25°C.

(Use a 30PA bottle (C) and 70PA bottle (C) at a temperature between 4°C and 20°C.)

#### SS; Stainless steel

Excellent in chemical resistance and heat resistance. It can be autoclaved at 121°C or less. For chemical resistance, see chemical resistance chart in separate manual.

## 2.3 Washing

Refer to table 2-1 for washing tubes and bottles.

- 1) Wash with tap water or dilute neutral detergent after use.
  - If there is a heavy stain, soak in dilute neutral detergent for a while or wash with a soft brush carefully.
- 2) Rinse with tap water and distilled water.
- 3) Dry in the air.

The PC tubes and bottles are inferior in chemical resistance with respect to alkali solution.

Do not use detergent of pH9 or above. And do not allow the soaking in dilute detergent for a long time because otherwise the material would become brittle.

For caps, take to pieces and wash in the same way with tubes. After washing, dry completely and store caps as a unit.

Take care of pH of detergent when performing ultrasonic washing.

Table 2-1 Washing Conditions

○: Usable, x: Not usable

	Tubes and	d bottles	Caps, Space caps,		Screw Caps of		
Washing conditions	PA,PP,PE,PET,SS	PC	Crowns, Spacers, Plugs	Adapters	bottles, O-ring		
Washing solution (pH5 or less)	×	×		×			
Washing solution	0	0					
(between pH5 and pH9)	O	0	0				
Washing solution	0	×					
(pH9 or above)	O	×	×				
Hot water (50°C or less)	0	0	0				
Ultrasonic washing	0	0		0			
Neutral detergent (pH7)	O	0	0				
Drying condition	Drying in the air						

⚠ CAUTION: In all cases, use neutral detergent (between pH5 and pH9).

#### 2.4 Sterilization

Refer to Table 2-2 for the sterilization of tubes and caps

#### How to autoclave tubes and bottles

- 1) Wash tubes and bottle well.
- 2) Stand it in the heat-resistance rack. Get off screw caps and inner caps of bottles to prevent deforming.
- 3) After autoclaving, take out tubes and . bottles after the temperature in the tank of autoclave lower to room temperature.

Table 2-2 Sterilization conditions

○: Satisfactory 🗙: Unsatisfactory

								Caps Space caps			Caps for bottles		
St	Sterilizing condition			Tubes, Bottles					rowr		Metal	Plastic	Adapter
		PA	PC	PET	PE	PP	SS	AL	TI	SS	caps	caps O-rings	
	115°C (0.7kg/cm²) for 30 min.	0	0	×	×	0	0	0	0	0	0	0*	0
Autoclaving	121°C (1.0kg/cm²) for 20 min.	0	0	×	×	0	0	0	0	0	0	0*	0
	126°C (1.4kg/cm²) for 15 min.		×	×	×	×	×	×	×	×	×	×	×
Boiling	15-30 min.	0	0	×	×	0	0	0	0	0	0	0	0
Ultraviolet rays	200-300 min	×	×	×	×	×	0	0	0	0	0	×	×
Gas	Ethylene oxide	×	×	×	0	0	0	0	0	0	0	0	0
Gas	Formaldehyde	0	0	0	0	0	×	0	0	0	0	0	0
	Ethanol (70%)	×	×	0	0	0	0	0	0	0	0	0	0
Chemical solution	Hydrogen peroxide (3%)	0	0	0	0	0	0	0	0	0	0	0	×
	Formalin	0	0	0	0	0	×	×	0	×	0	0	0


<sup>\*</sup> Green screw caps of bottles - old type - cannot be autoclaved. Blue one can be autoclaved.

## 2.5 Storage

Store the tubes and caps at room temperature. Do not store at high temperature or humidity or in a chemical vapor or expose to UV radiation.

## 2.6 Inspection

The cap, space cap and O-ring are consumable supplies. their lives vary in different conditions of use. Inspect them each time they are used, and replace them if you find any damage, deterioration or wear.



#### 2.7 Life

The life of plastic tubes and bottles depend upon the sample, speed of rotor, temperature, etc. The standard life of plastic tubes and bottles is specified as follows:

PET tubes, PA seal tubes, re-seal tubes, and cone top tubes are used one time only.

Table 2-3 Standard life of plastic tubes

			a o o. p.o.o o too o		
		Used at * standard condition	Used at over 70,000 rpm	Autoclaving	Used with weak** alkali sample
	PA, PP	5 times	1 times	5 times	5 times
Tubes	PC	5 times		5 times	5 times
	PE 5 times				5 times
Thick-walled tubes	,			50 times	50 times
Bottles	PC	20 times		5 times	10 times

<sup>\*</sup> In case they are used repeatedly for 24 hours at the maximum speed of rotors, below 70,000 rpm, with the aqueous sample.

Tubes capped with S series caps are disposable.

For caps, space caps, adapters and crowns, check for every use and replace worn or deteriorated one.

⚠ CAUTION: Do not use tubes/bottles that have exceeded their life expectancy. Failure to do so could result in damage of tubes/bottles and the rotor and the centrifuge.


The life expectancy of tubes/bottles depends on factors such as the characteristics of samples, speed of the rotor used, and temperature.


Always check for deterioration and damage (cracks, deformation, and so on) on tubes/bottles before using them. Do not use the tubes/ bottles if you find such a problem.

<sup>\*\*</sup> pH 7 to pH 9

## 2.8 How to use tube and bottles


• Tubes, 5mL and over, for angle rotors and vertical rotors.





# 2 Assembling the cap


- Assemble the cap loosely. Unscrew the setscrew.
- (2) Fill the tube upto three fourths with the sample, depress the cap until the tube comes in contact with the crown.
- (3) Tighten the nut manually.





# 3 Tighten the Cap

- (1) Insert the tube and the cap into the hole of the tube vice.
  - If there is a step between the top surface of the crown and that of tube vise, depress the cap until there is no step between them.
- (2) Tighten the handle (A) of the tube vise to fix the tube and cap. When you use the tube vise with the stopper, turn the cap to introduce the stopper into the slit of the stem on the cap.

If you tighten the nut at wrong position, the stopper may be damaged. Examples of correct position and wrong position on setting the stopper are shown in the below figures.







(3) Tighten the nut with a torque wrench or a box wrench.
If using the tube vise with the stopper, tighten the nut with a torque wrench with a clutch as holding down the stopper by your finger.

## ●Torque wrench

Torque wrench

## Torque wrench with a clutch



Tighten to function the clutch of the torque wrench.

# Socket

Tighten to a specified tightening torque.

#### Box wrench



Tighten the nut securely.


Tightening tools are specified for different caps (see appendix 4).

Table 2-4 Specified tightening torque

Part No. of cap	Name of cap (Volume of applied tube)	Tightening Torque
339673A	E-T12 Cap (40mL tube)	
S305231A	E3-AL Cap (40mL tube)	6 N⋅m
S308090A	E4-AL Cap (40mL tube)**	
322690A	F-AL Cap (94 mL tube)	0 NI
325968A	F2-AL Cap (94 mL tube)	8 N·m
463577A	C-TI Cap (12 mL tube)	
474070A	B2-TI Cap (6.5 mL tube)	10 N·m
481649A	C2-TI Cap (12 mL tube)	
474272A	E-TI Cap (40 mL tube)	12 N·m

\* 1N·m = 10.2kgf·cm


(4) Loose the handle (A) of the tube vice, and remove the tube with the tube setter.





<sup>\*\*</sup> E4-AL Cap can be used only for a 40PE tube.

# Balancing -


4








(2) Keep the difference of any two tubes which are arranged symmetrically within 0.3g. See 2.9. How to balance.








When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample.

Tighten the setscrew to prevent the leakage of the sample.

# •Tubes, 4mL or less, for angle rotors









## 4 Balancing

Put the tube and the cap in the balance, keep the difference of any two tubes which are arranged symmetrically with in 0.3g. When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample.




# 5 Fitting the cap

(1) Fit the cap to the tube depressing the cap until the tube comes in contact with the crown.




(2) Tighten the nut of the cap with the box wrench.



# 6 Tighten the Cap

- (1) Insert the tube and the cap into the hole of the tube vice and fix them.
- (2) Tighten the nut of the cap securely with the box wrench. If tightening is not securely, leakage will occur.
- (3) Remove the tube from the tube vice.




## ● Tubes (12PA/12PE/40PA/40PE tubes) that are applicable for S-Cap

S-Caps are applicable for P90AT, P80AT,P70AT2,P70AT, P50AT2 angle rotor. For more details, see the instruction manual (Part No. S999884) of S-caps.

⚠ CAUTION: If there is any abnormality such as corrosion, stop use of the S-Cap and contact a Hitachi Koki authorized sales/service representative.


NY Insert (Part No.:S401791A) are consumable parts. It is necessary to replace
it with the new one if it is worn. Contact a Hitachi Koki authorized sales/service
representative when replacing it.





# 1 Preparing the caps and adapters

Setscrew

① Confirm the composition of the cap. For about how to inspect the cap, refer to the instruction manual attached to the centrifuge or the rotor. ② Select the S-Cap, Tool Set, and Tube Holder according to the kinds of tubes.










Tube Holder C Tube Holder D (For 12 mL Tube) (For 30 mL Tube)

Tube Holder E (For 40 mL Tube)

Use the Tube Holder C inserted in Tube Holder E.








# 3 Depressing the Ti Ring


① Remove the tube with the Al Stem from the Tube Holder.

Then mount the Ti Ring on the tube from the bottom of the tube with the groove of the ring facing downward.




Ti Ring
The groove of the ring should be set downward.

② Set the tube as the above ① to the Tube Holder.




③ Place the Tube Holder on the suitable position of the Hand Press. Then push down the Handle until the Ti Ring contacts the step of the Al Stem.





# 4 Balancing


 Fill up the tube with the sample by using the injector etc.
 Otherwise the tube may be deformed during operation.



② Put the tube and setscrew in the Balance. Keep the difference of any two tubes which are arranged symmetrically within 0.3g.








<u>Tighten the setscrew to prevent the leakage of the sample.</u>

When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample.


#### Bottles



## filling the sample

Fill the bottle with the sample.

When a 30PC bottle (C3), 30PC bottle (C), 30PA bottle (C), 70PC bottle (B), 70PC bottle (C), 70PA bottle (C), 80PC bottle (C2), or 80PC bottle (C) are used at 100,000 x g or more (max. G-Force), the liquid should be over the shoulder of the bottle. Any volume which does not exceed the net volume is applicable at less than 100,000xg.



# 9 Balancing

Put the bottle and the cap in the balance, keep the difference of any two bottles which are arranged symmetrically within 0.3g. When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample within 5mm.




NOTE: Although you balanced the bottles within the approximate allowable imbalance (5mm), alarm message "IMBALANCE" might appear. Balance bottles more accurately again if the alarm message "IMBALANCE" appears (see the rotor instruction manual).

# 3 Tightening the cap


Fit the cap to the bottle and tighten the cap securely by hand.

#### Seal tubes



# filling the sample

Fill the samples of the same density into the tubes with an injector so a small quantity of bubbles remains. If the volume of the sample is small or the tubes are not welded completely, leakage will occur from the shoulder part or the tube may be deformed.




# 2 Balancing

Put the bottle and the cap in the balance, keep the difference of any two bottles which are arranged symmetrically within 0.3g. When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample within 5mm.




# 3 Sealing the tubes


- (1) Wipe off moisture adhering inside the opening of the seal tube.
- (2) Weld the seal tube according to the instruction manual of the STF-1, STF2 or STF3 tube sealer.




(i) Move the tube rack to the direction of arrows shown in the right figure to adjust position of the neck of the tubes (under the heater).




(ii) Depress the handle until the heater presses the neck of the tube. and then keep this condition for 2 to 3 seconds. Then depress the handle to the end and wait until the HEAT lamp is put out.



(3) Remove the protrusion on the welding part with your nail. Push the seal tube slightly and check that the sample does not leak. If the tubes are not welded completely, the tube may be deformed and the samples will leak.




## Thick-walled tubes for angle rotors



# filling the sample

Fill the tube with the sample.

You can use the thick-walled tube with optional volume less than net capacity.



# 2 Balancing


Put the tube in the balance. Keep the difference of any two tubes which are arranged symmetrically with in 0.3g.

When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample within 5mm.




↑ CAUTION: When using a 4.7PC thick-walled tube at the speed of over 60,000 rpm in the P100AT/P100AT2 rotor, cap a 4.7PC thick-walled tube with B-TI lid (Part No. S408523).

### Tubes for swing rotors



### filling the sample

Fill the tube up to within 3mm from the top of the tube with the sample.



### 2 Balancing


Put the tube in the balance, keep the differen6e of any two tubes which are arranged symmetrically within 0.2g.

When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample.

When using three tubes, keep the difference of three tubes within 0.2g.




### Re-seal tubes and cone-top tubes




### 1 Preparing the cap

- (1) Fix the tube vice to a desk, etc. by turning the handle (B).
- (2) Confirm the combination with the tube and the Crown assembly.



# Pilling the sample


Fill the samples of the same density into the tubes with an injector so a small quantity of bubbles remains. If the volume of the sample is small or the tubes are not welded completely, leakage will occur from the shoulder part or the tube may be deformed.



### 3 Balancing


Put the tube in the balance. Keep the difference of any two tubes which are arranged symmetrically with in 0.3g.

When using CP-NX, CP-WX, CP-MX, and CP- $\alpha/\beta$  series centrifuges, you can balance the tubes simply by leveling the sample within 5mm.



# 4 Attaching the crown and the plug to the tube

- (1) Insert the crown into the top of the tube.
  Push the crown to the top of the tube and <u>the bottom of the crown contact closely.</u>
- (2) Set the assembly of the tube and crown to the tube vice as shown in Fig. Be sure that there is no difference in level between the surface R of the crown and the surface (B) of the tube vice. If the surface (A) is higher than the surface (B), push the crown down to be no difference in level between (A) and (B).




- (3) Fix the crown and the tube by turning the handle (A) of the tube vice.
- (4) Tighten the plug from the top of the tube using a tube setter (C).

Normal assembly completion state is shown in Fig. (a).

To prevent the sample leakage during the operation it is important that the plug is securely fastened as shown in Fig.(a).

If the plug is not fastened securely as shown in Fig. (b), fasten it tightly.

Even though the plug is slipped by hard fastening, it can be used.




### HOW TO REMOVE THE TUBES FROM THE ROTOR

Draw out tubes from the rotor or the adapter with removal tool.

Before collecting samples, set the assembly of the tube and the crown to the tube vice again. And then remove the plug with the tube setter (C).

When removing the crown from the tube after collecting the samples, use the tube remover provided with the tube vice as shown in Figure.



#### 2.9 How to balance

Balance any two tubes with balance which sensitivity is less than 200mg, accessory of centrifuge. When using three tubes at swing rotor, Balance each of three tubes.

But using CP-NX, CP-WX, CP-MX, CP- $\alpha/\beta$  series centrifuge, you can balance the tubes simply by leveling the sample.

A CAUTION: Extremely different sample quantities must be avoided. Failure to do so might damage mechanical components.

⚠ CAUTION: Do not load only one tube or load tubes asymmetrically: Asymmetrically loading may cause imbalance operation and damage the centrifuge and rotor.




Fig. 2-6 How to balance

⚠ CAUTION: Be careful that imbalance operation may occur in the following cases.
 Fill the same sample in the tubes/bottles and load them in the rotor/buckets that are placed symmetrically with respect to the drive shaft in the rotor.
 If samples that are equal in volume but different in composition are used, the precipitation levels may be different by centrifugation and such operation may increase the level of imbalance.
 If samples that are equal in weight but different in volume (density) are used or if the tubes/bottles are different in inside diameter or shape, there may be variations in position of center of gravity and such operation may cause imbalance.

The tubes of the combination in Fig. 2-7 shall not set symmetrically because occurring imbalance during the operation.




Fig. 2-7 Correct balancing

### 3. How to use rotors

### 3.1 Allowable rotor speeds

The rotor should never be used at any speed higher than the maximum speed marked on the rotor in any 6ase. The rotor speed is limited by the following conditions:

DO NOT exceed this allowable speed.

#### (1) Allowable speed for density of the sample

To centrifuge a sample of an average density more than 1.2g / mL;
 Determine the allowable speed from the following equation.

Allowable speed (rpm) = maximum speed (rpm) 
$$\times \sqrt{\frac{1.2}{\text{Average density of a sample (g/mL)}}}$$

To centrifuge a sample with seal tubes in neo-angle rotors or vertical rotors;

These rotors can be used at their maximum speed with a sample with an average density less than 1.7g/mL. But in case of centrifuging a sample of an average density more than 1.7g/mL in these rotors, determine the allowable speed from the following equation.

Allowable speed (rpm) = maximum speed (rpm) 
$$\times \sqrt{\frac{1.7}{\text{Average density of a sample (g/mL)}}}$$

#### (2) Allowable speed for density gradient medium

The cesium chloride (CsCl) solution is frequently used as a density gradient medium, but if the CsCl solution with a high density is used, it may saturate during the rotation of the rotor depending on the rotation condition which may cause CsCl crystals to be educed. Educed crystals have a high density (approx. 4g / mL) and apply an excessive load to the rotor which is very dangerous. The crystallization will also greatly affect the density gradient of the CsCl solution and separation status of the sample. Therefore, always use a solution which will not crystallize within the allowable speed.

#### (3) Allowable speed for the combination of tubes, caps, etc..

The speed of the rotor should be limited by the combination of tubes, caps, etc.. For detail, consult "Rotors, Adapters, Tubes, Bottles and Caps", accessories of centrifuges.

### (4) Allowable speed for the life of aluminum rotors

The aluminum rotors should be limited the speed after 1,000 runs or 2,500 hours used See "4.6 Life of rotors."

### 3. 2 Sample limitation

> Your ultracentrifuge and rotor are not designed to confine any sample particles dispersed due to a leakage.


Therefore, when using radioactive, toxic or pathogenic materials, take additional precautions to prevent exposure to these materials, (e.g., use of isolated areas.)


### 3. 3 Adhesion of sample etc. .....

⚠ CAUTION : If sample etc. adheres to the rotor, use a soft cloth, etc. to wipe it off; leaving it could corrode the rotor.


Especially if sample adheres to the RLM adapter or optical adapter of rotor, immediately wipe it off; The adapter is aluminum alloy and susceptible to corrosion. If the RLM adapter or optical adapter corroded, the rotor will be unusable.


### 3.4 Angle rotors





- (1) Set the tube or the bottle in the rotor with the tube setter or the bottle setter.
- (2) When using the seal tube.
  Place the space cap on the tube to prevent deformation of the tube during operation.





CAUTION: Set the tube or the bottle symmetrically within the allowable imbalance, 0.3g or less.


# 3 Setting the cover


Fit the cover to the rotor and tighten the handle. Unless tightened adequately, packings may come off during operation.

A CAUTION: Be sure to attach the rotor cover to the rotor body and tighten handle securely.

Insufficient tightening of the handle can cause coming off the rotor cover and damage the centrifuge and the rotor.

• If using the rotor which handle separate into the handle and the handle bar, insert the handle bar to the hole of the handle.






### 4 Operation

(1) Gently and securely set the rotor onto the drive shaft.

⚠ CAUTION : If using the rotor which handle separate into the handle and the handle bar, be sure to remove the handle bar after setting the rotor.

(2) For operation, refer to the instruction manual of the preparative ultracentrifuge.

Drive shaft




# 5 Removal of sample

- (1) Upon completion of centrifuge, carefully remove the rotor from the drive shaft.
- (2) Loosen the handle and remove the cover.
- (3) Take tubes or bottles out with the tube setter etc.

  When using the seal tube under 10°C, it may be hard to take the tube out as soon as the rotor has stopped. In this case, take the tube out after 10 30 minutes.
- (4) Recover supernatants or pellets.
- - •Perform maintenance and inspection of the rotor each time it is used. If there is any abnormality, immediately stop using the rotor, and contact a Hitachi Koki authorized sales/service representative: Refer to "4.2 Maintenance of Rotors" and "4.4 Inspection of Rotors".


### 3.5 Swing rotors

Swing rotors (hook type buckets or pin type buckets) other than the P32ST rotor




### Preparation of the bucket

Remove the cap from the bucket.
 Stand the bucket in the bucket stand.




(2) Set the tube to the bucket. Be sure to set balanced tubes symmetrically.



(3) Tighten the cap of the bucket.

#### ●P65ST

Set the supplied driver into the slot of the cap and tighten it completely.



#### ●P28S, P28S2


Tighten the numbered cap to the same numbered bucket by hand to match each mark of the bucket and the cap.


Tighten the cap until the bottom of the cap bumps to the bucket and the cap must be adjusted to the match mark of the bucket.

#### ●P55ST2, P50S2, P40ST, P56ST

Tighten the numbered cap to the same numbered bucket with the supplied cap spanner to match each mark of the bucket and the cap.


Tighten the cap until the bottom of the cap bumps to the bucket and the cap must be adjusted to the match mark of the bucket.






### ? Setting the bucket

⚠ WARNING: When using swing rotor be sure to set all buckets whether or not samples are Put in buckets: Failure to do so could not only cause the rotor to vibrate, but could result in the rotor being deformed and the buckets becoming detached, which is very dangerous. Never use buckets made by other companies, or any other type of bucket that is not exclusively made for the rotor even if it is made by





Match numbered buckets with number of rotor and set the buckets to the rotor.



Hitachi Koki.

- (1) Insert the pin through holes of the rotor and the bucket.
- (2) Tighten them with supplied driver.






Hang the hook of the bucket on the pin of the rotor. When using the rotor which has two hooks, make sure that both hooks hang on the pin.

### ∠ Operation

★ WARNING: Before setting a rotor in ultracentrifuge, make sure that buckets have been securely installed in the rotor: Incorrect installation of buckets in rotor could damage the rotor or detach the buckets, which is very dangerous.


- •When using the rotor with hook-type buckets, be sure to hang on the bucket on the pin of the rotor correctly with a mirror or like this.
- Carry the large rotors with the supplied rotor setter.
  After setting the rotor, be sure to remove the rotor setter.
- (1) Gently and securely set the rotor onto the drive shaft.
- (2) For operation, refer to the instruction manual of the preparative ultracentrifuge.



# 5 Removal of sample


- (1) Carefully remove the rotor from the drive shaft.
- (2) Remove the buckets from the rotor.
  - ●P65ST · · · Unscrew the pin with the driver.
  - All rotors except P65ST · · · Remove the bucket by hand.
- (3) Remove the cap from the bucket with the driver or the cap spanner.
- (4) Take tubes out with pincettes.
- (5) Recover supernatants or pellets.
- CAUTION: Each time the rotor with optical adapter is used, be sure to enter the results in the "rotor log book": Failure to do so will disable management of the life of rotor. (No entry is needed for the rotor with RLM adapter.)
  - Perform maintenance and inspection of the rotor each time it is used.
     If there is any abnormality, immediately stop using the rotor, and contact a Hitachi Koki authorized sales/service representative: Refer to "4.2 Maintenance of Rotors" and "4.4 Inspection of Rotors"

#### ●P32ST rotor (top loading type buckets)




# Preparation of the bucket

Remove the cap from the bucket.
 Stand the bucket in the bucket stand.




(2) Set the tube to the bucket. Be sure to set balanced tubes symmetrically.



(3) Tighten the cap of the bucket.


Tighten the numbered cap to the same numbered bucket by hand to match each mark of the bucket and the cap.

Tighten the cap until the bottom of the cap bumps to the bucket and the cap must be adjusted to the match mark of the bucket.



### 3 Setting the bucket


⚠ WARNING: When using swing rotor be sure to set all buckets whether or not samples are Put in buckets: Failure to do so could not only cause the rotor to vibrate, but could result in the rotor being deformed and the buckets becoming detached, which is very dangerous. Never use buckets made by other companies, or any other type of bucket that is not exclusively made for the rotor even if it is made by Hitachi Koki.




Match numbered buckets with number of rotor and set the buckets to the rotor.

### 4 Operation

★ WARNING: Before setting a rotor in ultracentrifuge, make sure that buckets have been securely installed in the rotor: Incorrect installation of buckets in rotor could damage the rotor or detach the buckets, which is very dangerous.




- (2) After mounting the rotor on the drive shaft, check that the buckets swing smoothly by touching buckets gently.
- (3) For operation, refer to the instruction manual of the preparative ultracentrifuge.



# 5 Removal of sample


- (1) Carefully remove the rotor from the drive shaft.
- (2) Remove the buckets from the rotor.
- (3) Remove the cap from the bucket.
- (4) Take tubes out with tweezers.
- (5) Recover supernatants or pellets.
- ↑ CAUTION: Perform maintenance and inspection of the rotor each time it is used. If there is any abnormality, immediately stop using the rotor, and contact a Hitachi Koki authorized sales/service representative: Refer to "4.2 Maintenance of Rotors" and "4.4 Inspection of Rotors"

### 3.6 Neo angle rotors and titanium vertical rotors




### 2 Setting the tubes

(1) Load the rotor symmetrically with the tubes within allowable imbalance (0.3g or less) by using the tube setter, the removal tool or the tube setter (c). Excessive imbalance may damage the rotary shaft of the centrifuge.



(2) Place the space caps on the tubes to prevent deformation of the tubes during operation.




# $\mathcal{J}$ Setting the caps to rotor

- Using SEAL TUBE in P100VT, P83VT, P65VT2
- (1) Fix the holder base to a desk, etc..
- (2) Fix the rotor to the ditch of the holder base.
- (3) Place the cap on the rotor and tighten to 12N·m with the torque wrench.
- Using SEAL TUBE in P50VT2, P65VT3, P90NT, P65NT, P65NT2
- (1) Fix the projection of the handle to the hole of the cap.
- (2) Tighten the cap securely with the handle.
- (3) Remove the cap handle.

  Don't operate the rotor with the cap handle.
- Using RE-SEAL TUBE, CONE-TOP TUBE
- (1) Place the cap packing and the cap on the in tube in the rotor.
- (2) Tighten to 8N·m with the cap handle and the torque wrench.

The cap pushs down the tube is sealed. So tighten the cap securely to prevent leakage.

\* 1N·m= 10.2Kgf·cm



Torque wrench

\_ Cap

Cap

packing


2

Socket adapter

Hexagon socket

Rotor

↑ CAUTION: If the number of tubes to be set is small, do not install the space cap or the cap in the tube cavity into which no tube is inserted, but stay the cavities uninstalled for operation.




### 4

#### Operation

- (1) Gently and securely set rotor onto the drive shaft.
- (2) For operation, refer to the instruction manual of the preparative ultracentrifuge.

Longer acceleration and deceleration will result in better separation. If you start with homogeneous solution, slow deceleration only.

For more detail, refer to the instructions of the centrifuge.



# 5 Removal of sample

- (1) Upon completion of centrifuging, carefully remove the rotor from the drive shaft.
- (2) Remove the cap with the handle or the torque wrench.
- (3) Take space caps and tubes out with the tube setter.
- (4) Recover supernatant or pellets.

↑ CAUTION: • Each time the rotor with optical adapter is used, be sure to enter the results in the "rotor log book": Failure to do so will disable management of the life of rotor. (No entry is needed for the rotor with RLM adapter.)

Perform maintenance and inspection of the rotor each time it is used.
 If there is any abnormality, immediately stop using the rotor, and contact a Hitachi Koki authorized sales/service representative: Refer to "4.2 Maintenance of Rotors" and "4.4 Inspection of Rotors".

### 4. Maintenance

### 4.1 Anti-corrosion property

The anti-corrosion properties depend on the material of the rotor. See chemical resistance chart in separate manual.

- \* The over-speed adapters of all rotors are made from aluminum alloy.
- \* The covers of some titanium angle rotors are made from aluminum alloy. Consult Fig.1-1.
- \* The caps of all titanium vertical rotors are made from aluminum alloy.

The surface of aluminum rotors are anodized.

However, when salt solution (cesium chloride) is used in a swing rotor's bucket made from aluminum alloy, the bucket may be broken during centrifuging it it is corroded. For salt solution, use a rotor with a bucket made from titanium alloy or replace with a bucket made from titanium alloy. Consult the sales agent when replacing the bucket.

#### 4.2 Maintenance of rotors

Do maintenance the rotor to avoid corrosion after use.

#### (1) Ordinary maintenance

- 1) Wash the rotor with tap water or dilute neutral detergent.
- 2) Rinse with distilled water.
- 3) Wipe off water drops with soft cloth and dry the rotor up-side down.

After drying, coat the rotor with the silicone grease (vacuum grease). For drive hole, coat with silicone grease (vacuum grease) and wipe it. Store the rotor in a dry place.

Be sure to keep the RLM rotor on the rotor stand provided to protect the RLM adapter.




Fig.4-1 Coating the grease

For the rotor body of the swing rotor, wipe and coat with the silicone grease (vacuum grease).

#### (2) After using corrosive samples

Wash the rotor with flowing water for a while and then perform ordinary maintenance.

#### (3) If foreign matter adheres to the rotor

Soak the rotor in warm water (40~50°C) for a while, and remove foreign matter with a cleaning bar or a soft brush.

If you cannot remove it, contact a Hitachi Koki authorized sales/service representative.

MARNING: Do not allow the temperature of aluminum rotors and titanium rotors to rise above 100°C.

This would cause the material to become brittle.

- - Clean the inside of the drive hole (crown hole) of the rotor and the surface
    of the drive shaft (crown) of the centrifuge once a month. If the drive hole
    or the drive shaft is stained or any foreign matter is adhered, the rotor may
    be improperly installed and come of during operation.




Fig.4-2 How to the use a cleaning bar

### 4.3 Sterilizing rotors

Sterilize the rotor appropriately following Table 4-1.

MARNING: Do not autoclave the rotor or sterilize it in boiling water because its material may deteriorate and its strength may decrease.

Table 4-1 Sterilizing of rotor

			○:Usable	x:Not usable
		Aluminum rotors	Tita	anium rotors
	115°C (0.7kg/cm²) for 30 min.	×		×
Autoclaving	121°C (1.0kg/cm²) for 20 min.	×		×
	126°C (1.4kg/cm²) for 15 min.	×		×
Boiling	15-30 min.	×		×
Ultraviolet rays	200-300 nm.	0		0
0	Ethylene oxide	0		0
Gas	Formaldehyde	0		0
<b>a</b>	Ethanol (70%)	0		0
Chemical	Hydrogen peroxide (3%)	0		0
solution	3% Formalin*	0		0
He	Heat resistance (°C)			100


<sup>\*</sup>Do not dip the rotor in the formalin (3%) solution more than 2 hours.

## 4.4 Inspection of rotors


Check the rotor as follows after every use.

- · If the optical adapter/disk is unclean, wipe it with a soft cloth but do not damage it.
- · Replacements of packing and O-rings are available.
- · If the rotor appears to be corroded (deformation, cracks or discoloration is observed), do not use the rotor and call a Hitachi Koki authorized sales representative. The rotor may be irreparable depending on the degree of corrosion and wear.


### Angle rotors




#### Vertical rotors and Neo-angle rotors



#### Swing rotors




#### Buckets of swing rotors



↑ CAUTION: Check the packing (O-ring) every time. Replace it with new one if deterioration (crack, deformation, etc.) is observed. Otherwise, samples might be damaged due to poor rotor seal or the centrifuge or the rotor might be damaged during operation.

### 4.5 Replacing optical adapter/disk

The optical adapter/disk, if corroded or discolored, must be replaced immediately. Call the Sales Agent if it is necessary to replace the optical adapter. To replace the optical disk, follow the procedure below.



#### • Replacing the optical disk

#### (1) Preparation

- Prepare a new optical disk and knife.
   Make sure that the number of stripes of the new optical disk is right as meted in Table 4-2.
- · Wash and then dry the rotor well.




Table 4-2 Rotor's maximum speed and optical disk's number of stripes

Rotor's maximum speed (rpm)	Optical disk's number of black bands	Rotor's maximum speed (rpm)	Optical disk's number of black bands	Rotor's maximum speed (rpm)	Optical disk's number of black bands
100,000	9	50,000	18	30,000	29
90,000	10	48,000	19	28,000	32
80,000~83,000	11	45,000	20	27,000	33
70,000	13	42,000	21	25,000	35
65,000	14	40,000	22	23,000	39
60,000	15	38,000	24	21,000	40
55,000~56,000	16	35,000	25	19,000	45
54,000	17	32,000	28	17,000	50

- (2) Hold the rotor upside down.
- (3) Pry the edge of the optical disk with the knife and remove the disk. Be very careful not to damage the rotor.
- (4) Clean the disk hole on the rotor by using alcohol.
- (5) Remove the backing paper from the new optical disk.
  Place the optical disk so that it snugly fits into the groove of the disk hole on the rotor. Make sure the disk does not move.

#### • Pin Guide Replacement

When used with the 55PA, 65P and 55P-2 centrifuges, some rotors may need a pin guide (1.7 Relation between Rotor and Ultracentrifuge). If the rotors has a balance weight in place of a pin guide, replace it with a pin guide. If the pin protrudes from the pin guide, the pin guide should be replaced.




Fig.4-4 Pin guide attachment positions and how to attach

Use the spanner provided with the centrifuge to remove the pin guide or balance weight. Attach the new pin guide and tighten it with the spanner securely.

• When replacing the pin guide, check the number indicated on its head and use the guide the number of which corresponds to the maximum permissible speed of the rotor.

#### Example:

If the maximum speed of rotor is 50,000 rpm, the number should be 50. The maximum permissible speed decreases with an aluminum-alloy rotor that has exceeded its first lifetime. Be careful when selecting a pin guide for such a rotor.

### 4.6 Life of rotors

While using rotor repeatedly, its strength decreases gradually due to fatigue and creep of material, by which the rotor life is determined.

Swing rotors other than the P32ST rotor and Aluminum angle rotors

The swing rotors other than the P32ST rotor and aluminum angle rotors are assumed to be in the first lifetime until they reach 1,000 operation times or 2,500 accumulated hours. It is necessary to inspect rotors that have exceeded their first lifetimes at the works (this is charged).

Then the maximum permissible speed should be decreased by 10% and the rotor is assumed to be in its second lifetime until it reaches 1,000 operation times or 2,500 accumulated hours.

#### ●P32ST rotor

The P32ST rotor is assumed to be in the first lifetime until it reaches 500 operation times or 2,500 accumulated hours. It is necessary to inspect the rotor that has exceeded its first lifetime at the works (this is charged).

Then the maximum permissible speed should be decreased by 10% and the rotor is assumed to be in its second lifetime until it reaches 500 operation times or 2,500 accumulated hours.

The service life of RLM Rotor will be automatically calculated from the rotor speed and operating hours. Namely, if you use the rotor at a speed lower than its maximum permissible speed, the service life, in operating times and hours, will be prolonged automatically by calculating from the strength of rotor material.

●Titanium angle rotors, Neo-angle rotors and Vertical rotors

The titanium angle rotors, neo-angle rotors and vertical rotors must not be used when they reach 5,000 operation times or 10,000 accumulated hours. These rotors do not have a second lifetime.

For warranty, see the end of this manual.

■ When using the rotor with optical adapter, the "rotor log book" is essential to manage the life of rotor and for warranty when an accident occurs: Be sure to enter the results in it when using the rotor with optical adapter, and do not lose it. If any results of use are not entered in the "rotor log book", note that the rotor will not be covered by warranty.

# 5. Troubleshooting

#### 5.1 Rotors

Check the rotor immediately if it has something wrong.

Table 5-1 summarizes trouble symptoms, possible causes and actions to be taken.

• If the rotor has fallen on the floor by accident, have it checked by a Hitachi Koki authorized sales/service representative whether of not it is deformed.

If the rotor has fallen in the rotor compartment of the centrifuge, do not use the centrifuge and call a Hitachi Koki authorized sales/service representative immediately.

If there is only a trifle symptom of trouble, locate and remove the cause.
 If the cause cannot be located, do not use the device and call a Hitachi Koki authorized sales/service representative immediately.

Table 5-1 Troubleshooting of the rotor

Symptom	Possible causes and actions			
The rotor packing or cap packing breaks or comes out of position.	Are handle and cap tight securely?     Isn't packing elongated?     Is specified packing used?			
The handle or the cap cannot be tightened.	Is lubricant applied to threads?      →If threads are defective, call for repair.			
The rotor pin is bent.	→Call for repair.*			
The swing rotor's bucket hook is deformed.	Is the hook tightened at the wrong position set by the tightening tool?      →It should be repaired.			
The swing rotor's bucket does not swing normally.	<ul> <li>Is the bucket installed in reverse?(See 3.5 Swing rotors)</li> <li>Is the bucket cap tightened securely? Do the marks match?</li> <li>Is the bucket type correct?</li> <li>→If there is even a slight abnormality, inspect the rotor again.</li> </ul>			
The rotor is stuck on the shaft.	Dust or something like that is sticking to the shaft or drive hole.  Call a Hitachi Koki authorized sales/service representative.			
The aluminum rotor's tube cavity or bucket have changed color.	<ul> <li>Rinse sufficiently.</li> <li>If white dots are found after rinsing, inspection is necessary.</li> <li>Consult a sales agent.</li> <li>If there is too much corrosion, it may not be possible to repair it.</li> <li>If an aluminum bucket (old product) is used with a swing rotor, it could be corroded or broken unless it is maintained carefully. It is recommended to replace it with a titanium bucket.</li> </ul>			

Symptom	Possible causes and actions					
The ROTOR alarm lights on the centrifuge.	<ul> <li>The reflector (striped part) of the optical adapter or disk is unclean or damaged.</li> <li>→Wipe it with a soft cloth. Replace the adapter of disk if damaged.</li> <li>If the optical adapter needs replacement, call the Sales Agent. To replace the optical disk, see</li> <li>4.4 Inspection of Rotors ● replacing the optical disk</li> </ul>					
The IMBALANCE alarm lights on the centrifuge.	Is the sample well balanced?     Isn't the sample taking?					
The SPEED alarm lights on the centrifuge.	Isn't the speed set higher than the maximum permissible value?					
If another alarm lights on the centrifuge, refer to the operating instructions of the centrifuge.						

<sup>\*</sup> If the rotor needs to be repaired, it must be sent to the factory. So, call a Hitachi Koki authorized sales/service representative immediately.

Before shipping the rotor, be sure to sterilize and wash it.

#### 5.2 Tubes

The tube can be deformed or damaged by varying causes.

Check the following items in the order described. Additionally, see 2.6 Inspection and 2.7 Life ...

- (1) Isn't the amount of the sample too small?
  - $\rightarrow$  See 2.8 How to use Tube and Bottle

and 2.9 How to balance



- (2) Isn't the speed higher than the maximum permissible value?  $\rightarrow$

→ See 3.1 Allowable Rotor Speeds and the rotor instruction manual.



- (3) Are the cap and space cap, etc. combined properly?
  - Are the parts of the cap combined properly?
  - $\rightarrow$  See appendix 1, appendix 2, appendix 3, and the rotor instruction manual.









- (4) Is the cover or the cap of the rotor tight?
  - · Is the cap of the tube tightened firmly?
  - When using seal tubes, is the tube sealer adjusted properly?
     (heater height, centering) → See the STF-1 instructions, the STF2 instructions, or the STF3 instructions.

Are burrs removed after welding?  $\rightarrow$  See the STF-1 instructions, the STF2 instructions, or the STF3 instructions.

- · Is the tightening torque of the caps of the neo-angle rotors or the vertical rotors right?
  - → See 3.6 Neo angle rotors and titanium vertical rotors.





(6) Is the sample and detergent good?  $\rightarrow$  See 2.2 Materials of Tubes . Are sterilizing conditions right?  $\rightarrow$  See 2.4 Sterilization .

- If the tube is squeezed, fill the tube hole on the rotor with tepid water and, after a few minutes, pull out the tube by using tweezers or something like that.
- If there is only a trifle symptom of trouble, locate the cause or call a Hitachi Koki authorized sales/service representative.


If such a symptom is overlooked, an accident might result.


- If the cap has fallen in the bottom of the tube hole and you cannot removed it, call a Hitachi Koki authorized sales/service representative.
- \* When calling a Hitachi Koki authorized sales/service representative, keep the part of interest available for examination.

# Appendix 1 The list of caps for tubes


Different caps use a different composition or parts materials. Be careful not to combine unsuitable parts.

### Component





### Materials



	1					1			ī
Parts Caps	Set screw	Nut	Washer	Crown	Weight ring	O-ring	Gasket	Stem	Insert
A2-AL Cap 336710A		80130235 AL		474117 AL		S406640A NBR		474116 AL	
A3-AL Cap 336711A		80130058 AL		474123 AL		S401807A NBR		474121 AL	
B-AL Cap	S401829A SS	441393		414432			S401778A	453411A	S401791A
414429A B2-AL Cap	S401829A	AL 441393		AL 474641		S401821A	<u>CR</u>	AL 474639A	NY S401791A
474507A B-TI Cap	SS S401829A	AL 459561		AL 459559		NBR	S406626A	AL 460697A	NY S401791A
460696A B2-TI Cap	SS S401829A	TI 4595612		TI 474073		 S401815A	NBR	TI 474071A	NY S401791A
474070A	SS	TI		TI		NBR		TI	NY
C-AL Cap 441399A	S401829A SS	441393 AL		441392 AL			S401796A CR	453412A AL	S401791A NY
C4-AL Cap S407404A	S401829A SS	441393 AL		S407407 AL		S401803A NBR		S407405A AL	S401791A NY
C-TI Cap 463577A	S401829A SS	463567 TI		463569 TI		S401803A NBR		463576A TI	S401791A NY
C2-TI Cap 481649A	S401829A SS	4595612 TI		481652 TI		S401803A NBR		481650A TI	S401791A NY
C-SS Cap 413312A	S401829A SS	413309		413308 SS			S401796A CR	453413A	S401791A
E-AL Cap	S401829A	SS 434022	S401813A	435519	459814	S401787A	CR	SS 453416A	NY S401791A
321517A	SS	AL	POM	AL	AL	NBR		AL	NY
E2-AL Cap 331198A	S401829A SS	434022 AL		434021 AL	459814 AL	S401787A NBR		453416A AL	S401791A NY
E3-AL Cap S305231A	S401829A SS	S408401 AL	S401813A POM	S408400 AL	459814 AL	S401787A NBR		S408421A AL	S401791A NY
E4-AL Cap	S401829A	S408401	S401813A	S411239A	459814	S401787A	_	S408421A	S401791A
S308090A E-TI Cap	SS S401829A	AL 474277	POM	AL 474276	AL 474275	NBR S401787A		AL 474273A	NY S401791A
474272A	SS	TI		TI	TI	NBR	<b></b>	TI	NY
E-Tl2 Cap 339673A	S401829A SS	434421 AL	S401811A POM	478423 TI	459814 AL	S401787A NBR		478421A AL	S401791A NY
E-SS Cap 451480A	S401829A SS	4187622 SS		451481 SS			S406625A CR	453417A SS	S401791A NY
E2-SS Cap 418763A	S401829A SS	4187622 SS		418765 SS			S406625A CR	453417A SS	S401791A NY
F-AL Cap	S401829A	434421	S401811A	436127	459818	S401792A		453414A	S401791A
322690A F2-AL Cap	SS S401829A	AL 434421	POM S401811A	AL 435766	AL 459818	NBR S401792A		453414A	NY S401791A
325968A F-SS Cap	SS S401829A	AL 418762	POM	AL 428761	AL	NBR	S401797A	AL 453415A	NY S401791A
418759A	SS	SS		SS			CR	SS	NY

NOTE1) the upper = part No. the lower = materials


NOTE2) Set screw, washers, O-rings, gaskets and inserts are available in sets of 10.

NOTE3) Individual nuts, crowns, weights and stems are available. The stem is provided with an insert.

# Appendix 2 The list of S-Series Cap and tools

## S-Series Cap (For 12PA/12PE/40PA/40PE tubes)

### Component




Parts Caps	Set screw	AL Stem	Insert	Ti Ring
S-12AL Cap	S401829A	S410543A	S401791A	S410545
S410542A	0 <del>4</del> 01023A	AL	NY	TI
S-40ALCap	C404920A	S410533A	S401791A	S411473
S410532A	S401829A	AL	NY	TI

NOTE1) the upper = part No. the lower = materials

NOTE2) Setscrew and inserts are available in sets of 10.

NOTE3) Individual AL Stems and Ti Ring are available.

## S-Cap Tool Set



Cap		Applied rotors			
Сар	Tool set	Tube Holder	Tube Setter	Hand Press	Applied folois
S-12ALCap S410542A	•		481056A	S204591A	P90AT, P80AT,
01100121		S411486			170/112
S-40ALCap S410532A	• •		Tube Setter(C)	Hand Press	P70AT, P50AT2

\*Use the Tube Holder C inserted in Tube Holder E.

# Cap Tool Kit

Douto	Сар		Cap Tube holder E		Tube holder C		Tube setter (C)		Hand Press	
Parts  Kits	9									
	Part No./Name	Qty.	Part No./Name	Qty.	Part No./Name	Qty.	Part No./Name	Qty.	Part No./Name	Qty.
S308626A S-12AL Cap Tool Kit	S410542A S-12AL Cap	8	S411486 Tube holder E	1	S411487 Tube holder C	1	481056A Tube setter (C)	1	S204591A Hand Press	1
S308626B S-12AL Cap Tool Kit	S410542A S-12AL Cap	12	S411486 Tube holder E	1	S411487 Tube holder C	1	481056A Tube setter (C)	1	S204591A Hand Press	1
S308627A S-40AL Cap Tool Kit	S410532A S-40AL Cap	8	S411486 Tube holder E	1			481056A Tube setter (C)	1	S204591A Hand Press	1
S308627B S-40AL Cap Tool Kit	S410532A S-40AL Cap	12	S411486 Tube holder E	1			481056A Tube setter (C)	1	S204591A Hand Press	1

# **Appendix 3** The list of space caps


Select the correct space caps to match the tubes and rotors to be used.

	S	oace caps	Applicable rotors (Old models)		
Applicable tubes	Part No.	Name (color)			
2PA Seal tube	S402631	Space cap (A) (blue)	(RP70VT)		
	484284	B Space cap (white)	(RPV65T)		
5PA Seal tube	488101	B2 Space cap (black)	P100VT, P65VT2, P90NT, P65NT2, (P83VT), (P65VF), (P55VF2), (RP65AF), (SRP83VT), (RP65VF), (RP55VF2), (RP65VT2), (RP67VF), (RP85NT)		
	S403331	B3 Space cap (white)	(RP50AT4), (SRP50AT)		
6.5PA Seal tube	S407828	B5 Space cap (red)	P100AT2, P50AT4, (P100AT), (RP50AT4), (RP55), (SRP50AT), (RP40-3), (RP40-2)		
	484285	C Space cap (white)	P65A, (P65AT), (RP65T), (RP65), (RP40), (RPW65T), (RPW45), (P55AT)		
12PA Seal tube	486582	C2 Space cap (red)	P90AT, P80AT, P70AT2, (RP83T), (RP80T-3), (SRP70AT)		
	486583	C3 Space cap (white)	(RPV45T)		
	S402310	C4 Space cap (black)	P65VT3, P65NT, (RP55VF), (RP65VT3), (RP65NT)		
40PA Seal tube	485649	E Space cap (red)	P70AT, P50AT2, (P50A2), (P30A2), (RP70T), (RP50T-2), (RP50T), (RP50-2), (RP30-2), (RPW50T), (RP60T)		
	484252	E2 Space cap (blue)	(RPV50T), (RPV30)		
	S402525	E3 Space cap (white)	P50VT2, (RP50VF)		
94PA Seal tube	S408969	F-Ti Space cap (red)	P45AT, (P42A), (RP45T), (RP42), (RPW35)		


The tube sealer, model STF-1, model STF2, or model STF3, and the tube rack are necessary to seal the tube.

Applicable seal tubes	Tube racks	Part No.	Remarks		
1.5PA seal tube	Tube rack (G2)	S201778G			
2PA seal tube	Tube rack (G)	S201778F			
3.5PA seal tube	Tube rack (B2)	S201778E			
4PA seal tube	Tube rack (B3)	S201778H			
5PA seal tube	Tube rack (B)	S201778A	For STF-1/STF2/STF3		
6.5PA seal tube	Tube rack (B4)	S201778J			
8PA seal tube	Tube rack (C2)	S201778L			
12PA seal tube	Tube rack (C)	S201778B			
40PA seal tube	Tube rack (E)	S201778C			
O4DA applitude	Tube rack (F)	S201778D	For STF-1		
94PA seal tube	Tube rack (F2)	S201778M	For STF2/STF3		
2PA,4PA,5PA,6.5PA,12PA			For STF-1/ STF2/STF3		
and 40PA seal tubes	Multi-rack	S201778K	The 1.5PA, 3.5PA, 8PA and 94PA seal tubes		
			cannot be used in the multi-rack.		

Type of tube rack is stamped on top of the tube rack. e.g.) In case of Tube rack (G), "G" is stamped on top of the tube rack.



# Appendix 4 Tightening tools of the cap







<Box wrench>



<Torque wrench> #469283 #474871



<Torque wrench> #477456A

Caps				
(tightening	Tool sets	Parts of	the tool set	Applied rotors
torque kg/cm)	100/3613	Tube vises	Box wrenches Torque wrenches	(Old models)
A2-AL Cap		S303696A	S406641	
336710A		Tube vise A2	Box wrench	
A3-AL Cap		S303696A	S406641	
336711A		Tube vise A2	Box wrench	
	S300715A	348253A	4361453	P50AT4, (RP50AT4),
	Tool set (B)	Tube vise (B)	Box wrench	(SRP50AT)
B-AL Cap	, ,	336662C	4361453	·· · ·
414429A		Tube vise (B)	Box wrench	(For 5PET tube)
	S300716A	S300717A	4361453	(DD00.0)
	Tool set (B2)	Tube vise (B2)	Box wrench	(RP30-3)
B2-AL Cap	` ,	336662C	4361453	P50AT4, (RP50AT4)
474507A		Tube vise (B)	Box wrench	(For 6.5PET tube)
B-TI Cap	S300715A	348253A	4361453	
460696A	Tool set (B)	Tube vise (B)	Box wrench	P100AT2, (P100AT), (RP80T)
B2-TI Cap		336662A	469283	(DD) (CET)
474070A (100)		Tube vise (B)	Torque wrench	(RPV65T)
C-AL Cap	S300713A	S300714A	4361453	P65A, (P65AT), (P55AT), (RP65T),
441399A	Tool set (C)	Tube vise (C)	Box wrench	(RP65)
C4-AL Cap		333885B	4361453	P65A, (P65AT), (P55AT),
S407404A		Tube vise (C2)	Box wrench	(For 12PET tube)
C-TI Cap	341109A	333885A	469283	P90AT, P80AT, P70AT2,
463577A(100)	Tool set (C2)	Tube vise (C2)	Torque wrench	(SRP70AT), (RP83T), (RP80T-3)
C2-TI Cap		333885B	469283	(RPV45T)
481649A (100)		Tube vise (C2)	Torque wrench	,
C-SS Cap	S300713A	S300714A	4361453	P65A, (P65AT), (P55AT), (RP65T),
413312A	Tool set (C)	Tube vise (C)	Box wrench	(RP65)
E4-AL Cap	339198A	343417A	477456A	P70AT, P50AT2
S308090A (60)	Tool set (E)	Tube vise (E)	Torque wrench	(For 40PE tube)
E3-AL Cap	339198A	343417A	477456A	P70AT, (RP70AT), P50AT2,
S305231A (60)	Tool set (E)	Tube vise (E)	Torque wrench	(RP50T-2), (RP50T), (RPW50T),
` ′	` '	` '	•	(RP50-2), (P50A2), (P30A2), (RP30-2)
E2-AL Cap	339198A	343417A	477456A	(RP30)
331198A	Tool set (E)	Tube vise (E)	Torque wrench 474871	, ,
E-TI Cap		338956A		(P50VT), (RPV50T), (RPV30)
474272A (120)	339198A	Tube vise 343417A	Torque wrench	, , , , , , , , , , , , , , , , , , , ,
E-Tl2 Cap 339673A (60)	339198A Tool set (E)	343417A Tube vise (E)	477456A Torque wrench	P70AT, (RP70T)
	1001 SEL (E)	343417A	4361453	
E-SS Cap 451480A		Tube vise (E)	Box wrench	P70AT, P50AT2, (P50A2), (P30A2)
E2-SS Cap		343417A	4361453	
418763A		Tube vise (E)	Box wrench	
F-AL Cap	339083A	215786B	477456A	
322690A (80)	Tool set (F)	Tube vise (F)	Torque wrench	(RP21)
F2-AL Cap	339082A	215786A	477456A	
325968A (80)	Tool set (F2)	Tube vise (F)	Torque wrench	P45AT, (P42A)
F-SS Cap	100.001 (1.2)	215786B	4361453	(222.1)
418759A		Tube vise (F)	Box wrench	(RP21)
410/09A		Tube vise (F)	DOX WIEITCH	<u> </u>

# **Appendix 5** The list of bottles

Polypropylene copolymer


How to indicate Parts of bottles, Part No., and material

- The upper number is part No. The parts with ★ mark are 1 pc/set, and others are 10 pcs/set. If blank, the part is not provided.
- The lower letters indicate materials.

Material

PA: polypropylene copolymer PPO: polyphenylene oxide PP: polypropylene AL: aluminum alloy PC: polycarbonate POM: polyacetal

NBR: nitrile-butadiene rubber



Vol. * <bottles></bottles>		Pa	rts of bottles, Pa	Applied rotors		
(mL)	Name	Screw cap	Inner cap	O-ring	Bottle	(old models)
	334105A (10pcs/set) 10PA Bottle	S401802A PP			PA	P90AT, P80AT, P70AT2, (P65AT), (P55AT), (SRP70AT), (RP65T)
10	325952A (10pcs/set) 10PC Bottle	S401802A PP			PC	P90AT, P80AT, P70AT2, (P65AT), (P55AT), (SRP70AT), (RP65T)
	S303922A (6pcs/set) 10PC Bottle (B)	474042A PPO	474043A PPO	474044A NBR	336626A PC	P90AT, P80AT, P70AT2, (P65AT), (P55AT), (SRP70AT), (RP65T)
	330345A (10pcs/set) 30PA Bottle	S401789A PP			PA	P70AT, P50AT2, (P50A2), (RP70T), (RP50T-2), (RP50-2)
30	S308133A (6pcs/set) 30PC Bottle (B3)	474035A PPO	474036A PPO	474037A NBR	S308134A PC	P70AT, P50AT2, (P50A2), (RP70T), (RP50T-2), (RP50-2)
	S308132A (6pcs/set) 30PC Bottle (C3)	485541★ AL	474036A PPO	474037A NBR	S308134A PC	P70AT, P50AT2, (P50A2), (RP70T), (RP50T-2), (RP50-2)
	S310626A (6pcs/set) 30PA Bottle (C)	485541★ AL	S413364A PP	474037A NBR	S310625A PA	P70AT, P50AT2
	S308433A (6pcs/set) 70PC Bottle (B)	474039A PPO	474040A PPO	474041A NBR	S308431A PC	P45AT, (P42A), (RP45T), (RP42), (RPW35)
70	S308405A (6pcs/set) 70PC Bottle (C)	486370★ AL	474040A PPO	474041A NBR	S308431A PC	P45AT
	S310331A (6pcs/set) 70PA Bottle (C)	486370★ AL	S413195A PP	474041A NBR	S310419A PA	P45AT
230	S413248B (2pcs/set) AL Cap (2) assy **	S310455★ AL	474555★ POM	474633★ NBR		P21A2
230	S403617A (1pc/set) AL Cap assy **	S301384 <b>★</b> AL	S401808A POM	S401809A NBR		P19A, (RP19)

<sup>\*</sup> Nominal capacity


<sup>\*\*</sup> The caps for 230 PA bottles.

# Appendix 6 The list of tube kits

### ●Open-top Tube Kit

Parts	Parts TUBE		САР		PACKING for CAPS		TUBE VISI		WRENCH 4361453 469283	) )	TUBE SETTER		
	Part No./Name	**	Part No./Name	**	Part No./Name	**	Part No./Name	**	Part No./Name	**	Part No./Name	**	
S304985A 6.5PA Tube Kit	329445A 6.5PA Tube	2 (100)	460696A B-TI Cap	8	S406626A	2 (20)	348253A Tube Vise (B)	1	4361453 Box wrench	1	435823A Tube Setter	1	
S304986A 12PA Tube Kit (A)	329606A 12PA Tube	2 (100)	463577A C-Tl Cap	12	S401803A	2 (20)	333885A Tube Vise (C2)	1	469283 Torque wrench	1	435823A Tube Setter	1	
S304986B 12PA Tube Kit (B)	329606A 12PA Tube	2 (100)	441399A C-AL Cap	12	S401796A	2 (20)	S300714A Tube Vise (C)	1	4361453 Box wrench	1	435823A Tube Setter	1	
S304987A 40PA Tube Kit (A)	329607A 40PA Tube	2 (100)	339673A E-Tl2 Cap	8	S401787A	2 (20)	343417A Tube Vise (E)	1	*1	1set	435823A Tube Setter	1	
S304987B 40PA Tube Kit (B)	329607A 40PA Tube	2 (100)	S305231A E3-AL Cap	12	S401787A	2 (20)	343417A Tube Vise (E)	1	*2	1set	435823A Tube Setter	1	
S304988A 94PA Tube Kit	S304299A 94PA Tube	4 (100)	325968A F2-AL Cap	6	S401792A	1 (10)	215786A Tube Vise (F)	1	*1	1set	435823A Tube Setter	1	

<sup>\*\*</sup> Q'ty (pcs/box)



#### •Seal Tube Kit

Parts SEALTU		<b>=</b>	SPACE CAI	P	TUBE SETTER		
Kits							
	Part No./Name	**	Part No./Name	**	Part No./Name	**	
S304989A	345319A	2	488101	8	S407157	1	
5PA Seal Tube Kit (A)	5PA Seal Tube	(100)	B2-Spase Cap	0	Tube Setter	ı	
S304989A	345319A	2	488101	18	S407157	1	
5PA Seal Tube Kit (B)	5PA Seal Tube	(100)	B2- Spase Cap	10	Tube Setter	1	
S304990A	S304238A	2	S407828	8	S407157	1	
6.5PA Seal Tube Kit	6.5PA Seal Tube	(100)	B5- Spase Cap	0	Tube Setter	'	
S304991A	345320A	2	486582	12	S407157A	1	
12PA Seal Tube Kit (A)	12PA Seal Tube	(100)	C2- Spase Cap	12	Tube Setter	'	
S304991B	345320A	2	484285	12	S407157	1	
12PA Seal Tube Kit (B)	12PA Seal Tube		C-Spase Cap	12	Tube Setter	'	
S304991C	345320A	2	S402310	10	S407157	1	
12PA Seal Tube Kit (C)		`			Tube Setter		
S304992A	345321A	2	485649	12	S407157	1	
40PA Seal Tube Kit (A)				'-	Tube Setter		
S304992B	345321A	2	S402525	8	S407157	1	
40PA Seal Tube Kit (B)		` '			Tube Setter		
S304993A	S304310A	4	485650	6	S407157	1	
94PA Seal Tube Kit	94PA Seal Tube	(100)	F2- Spase Cap		Tube Setter	•	

<sup>\*\*</sup> Q'ty (pcs/box)

The tube sealer, model STF-1, or model STF2, model STF3 and the tube rack are separately necessary for using seal tubes. See appendix 3 for the part No. of the tube rack.

# Applicable rotors

Applicable rotors of every tube kit are shown below.

Applicable rotors of every tube kit are shown below.						•	• : Applicable										
Kits	Rotors	P100AT2	P90AT	P80AT	P70AT2	P70AT	P65A	P50AT2	P45AT	P32AT	P90NT	P65NT	P65NT2	P100VT	P65VT3	P65VT2	P50VT2
S304985A	6.5 PA Tube Kit	•															
S304986A	12 PA Tube Kit (A)		•	lacksquare	•												
S304986B	12 PA Tube Kit (B)						•			•							
S304987A	40 PA Tube Kit (A)					•											
S304987B	40 PA Tube Kit (B)							•									
S304988A	94 PA Tube Kit								•								
S304989A	5 PA Seal Tube Kit (A)										•			•			
S304989B	5 PA Seal Tube Kit (B)												•			•	
S304990A	6.5 PA Seal Tube Kit	•															
S304991A	12 PA Seal Tube Kit (A)		•	•	•												
S304991B	12 PA Seal Tube Kit (B)						•										
S304991C	12 PA Seal Tube Kit (C)											•			•		
S304992A	40 PA Seal Tube Kit (A)					•		•									
S304992B	40 PA Seal Tube Kit (B)																•
S304993A	94 PA Seal Tube Kit								•								

### Warranty

#### Warranty for Ultracentrifuge Rotors

Hitachi Koki warrants its ultracentrifuge rotors under the terms and conditions set forth below.

#### 1. Scope of Application

This warranty applies to:

- (1) failure or damage to the rotor due to rotor defects, in materials or workmanship, that are attributable to Hitachi Koki, and
- (2) failure or damage to a Hitachi Koki ultracentrifuge due to defects in a rotor covered by the warranty.

Note: This warranty covers the rotor or ultracentrifuge only, and Hitachi Koki shall not be liable for damage to samples or other items that may arise from failure of or damage to the rotor or ultracentrifuge.

#### 2. Warranty Period

This warranty is valid for five years from the date of delivery of the rotor by Hitachi Koki or a Hitachi Koki dealer. However, except as provided under (1) ②below, even within such period the warranty does not cover cases where the user has exceeded the applicable limit on the number of runs or the cumulative operating time for the type of rotor, as set forth in (1) and (2) below.

- (1) Aluminum alloy angle rotor, swing rotor, zonal rotor, and continuous flow rotor
  - ① The limits on the maximum number of runs and the maximum cumulative operating time for these rotor types are stated in the table below. The warranty covers the time up until either of these limits is reached, provided that it is reached during the warranty period and that the maximum and allowable speeds are not exceeded.
  - ② In order for the warranty to cover use after either of the limits in the table below is reached, inspection by Hitachi Koki will be necessary. The user must pay for such inspection. If Hitachi Koki's inspection finds that the rotor is free of corrosion and other abnormality, Hitachi Koki will then further warrant the rotor up until either of the limits in the table below is reached once more, but with the rotor's maximum speed reduced by 10%.

Even in such case however, the warranty period shall remain limited to five years from the date of delivery and shall not be extended.

	Maximum number of runs	Maximum cumulative operating time
P32ST rotor	500	2,500 hours
Rotors other than the P32ST rotor	1,000	2,500 hours

(2) Titanium alloy angle rotor, titanium alloy vertical rotor, and neo-angle rotor

The limits on the maximum number of runs and the maximum cumulative operating time for these rotor types are stated in the table below. The warranty covers the time up until either of these limits is reached, provided that it is reached during the warranty period and that the maximum and allowable speeds are not exceeded.

(With these types, the warranty cannot be extended beyond the stated limits, even if the rotor undergoes inspection by Hitachi Koki.)

Maximum number of runs	Maximum cumulative operating time
5,000	10,000 hours

#### 3. Warranty Terms

(1) Should a rotor covered by this warranty fail or be damaged, Hitachi Koki will replace it with a new rotor in return for payment of a sum determined according to the formula below, based on the ratio of the actual number of runs or cumulative operating time to the maximum number of runs or cumulative operating time.

(2) Should an ultracentrifuge fail or be damaged due to failure of or damage to a rotor covered by this warranty, Hitachi Koki will supply, free of charge, all parts required for repair of the ultracentrifuge. However, in the case of failure of or damage to an ultracentrifuge's drive unit, Hitachi Koki will replace it with a new drive unit in return for payment of a sum derived by multiplying the current price of the drive unit by the ratio of the actual total number of revolutions or years completed to the warranted total number of revolutions or warranted years that are stated in the operator's manual for the ultracentrifuge, as in the formula below.

#### 4. Notes

- (1) This warranty extends only to the original buyer from Hitachi Koki or the original buyer from a Hitachi Koki dealer. Rotors whose ownership passes from the original buyer to a third person due to resale, transfer or the like are not covered by this warranty. Neither does this warranty cover cases where the ownership of a rotor is reacquired subsequently by the original buyer.
- (2) The warranty for the rotor proper (warranty period and terms) shall continue to apply in cases where rotor parts and the like are additionally purchased. In such cases, the warranty period shall continue to be measured from the date of delivery of the rotor. If such parts are subject to a restricted warranty period, this warranty shall only cover failure or damage involving such parts that occurs both within the rotor warranty period and within the warranty period for the parts concerned.
- (3) This warranty does not cover any of the following, even during the warranty period:
  - ① Failure or damage due to abuse or misuse.
  - ② Failure or damage due to operation or maintenance in a manner contrary to the instructions in the operator's manuals for the rotor and ultracentrifuge.
  - 3 Failure or damage due to disassembly or modification without the permission of Hitachi Koki.
  - Failure or damage due to use of a bucket, adapter, tube, bottle, or the like, other than those
     approved by Hitachi Koki for use with the product.
  - ⑤ Failure or damage due to disaster such as fire or earthquake.
  - © Failure or damage due to use with an ultracentrifuge not of Hitachi Koki manufacture without the permission of Hitachi Koki.
  - Pailure or damage due to use with an ultracentrifuge that has been modified without the permission of Hitachi Koki.
  - Failure or damage due to a failed or damaged part (including a tube, bottle, adapter, etc) not submitted to Hitachi Koki.
  - Failure or damage not recorded in a rotor use log book (in the case of an optical adapter type rotor).
  - ® Replacement of consumable items.
  - 1 Parts subject to a restricted warranty period that has expired.

#### Other Matters

Rotors that are delivered to Hitachi Koki for inspection or repair must be sterilized and cleaned beforehand.

Be aware that Hitachi Koki will sometimes decline to carry out inspection or repair if a rotor or other related equipment is or could be contaminated as a result of use without protection or under inadequate safety management in a radioisotope facility or a laboratory of the level P2 or higher.

### **Decontamination**

#### **⚠** WARNING:

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from a Hitachi Koki authorized sales/service representative. Note that Hitachi Koki cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- ●It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to a Hitachi Koki authorized sales/service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. Hitachi Koki may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by Hitachi Koki. It is your responsibility to bear the cost of sterilization or decontamination. If you have any question, please send e-mail to "himac@hitachi-koki.co.jp". Note that Hitachi Koki cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

### Rotor retirement

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years		
Angle rotor	Titanium alloy	12		
Vertical rotor	Aluminum alloy	10		
Swing rotor	Titanium alloy Aluminum alloy	10		
Zonal rotor	Titanium alloy	10		
Continuous flow rotor	Titanium alloy	10		

# 

1060,Takeda,Hitachinaka City Ibaraki Pref., 312-8502 Japan

Tel: (81)29-276-7914 (Dial in)

Fax: (81)29-276-7475

URL: http://www.hitachi-koki.com/himac/